

# Life Scope *TR*

Bedside Monitors BSM-6000 series

## See the difference



*Fighting Disease with Electronics*

 **NIHON KOHDEN**

# Valuable features without additional cost


## Smart Cable™ system - new modular technology



Smart Cable™ technology miniaturizes circuits found in traditional modules and embeds that circuitry into the cable. With Smart Cable™ technology, you get complete modular flexibility at a significantly reduced cost and without the inconvenience associated with traditional modular systems.

When you plug a Smart Cable™ into a MULTICONNECTOR, it automatically detects the type of parameter and starts measuring.

As new parameters are developed, Nihon Kohden will incorporate the new parameters into the new Smart Cable™, guaranteeing your monitoring investment against obsolescence.



## Smart data review

Life Scope TR provides storage and review capabilities within the bedside monitor that are typically found only in a central station. Moreover review screens can be synchronized with each other. This enables you to easily access any necessary information.

### Full disclosure waveforms



Up to 72 hours of 5 selected waveforms.

### Arrhythmia recall



Over 8,000 arrhythmia events. Multi-template arrhythmia analysis provides more reliable arrhythmia information.

### ST template



Multi-lead ST segment monitoring provides you with continuous oversight of transient changes in your patient's cardiac condition.

### 12-lead ECG analysis




Life Scope TR includes ECAPS 12C, the same 12-lead ECG interpretation software in Nihon Kohden electrocardiographs.

# Intelligible Fluid Management

## esCCO - Non-invasive continuous cardiac output from ECG and SpO<sub>2</sub>

Nihon Kohden is redefining Quality of Care with new non-invasive technologies like PWTT and esCCO by introducing volumetric information to all care levels.


esCCO (estimated continuous cardiac output) provides real-time, continuous non-invasive cardiac output measurement alongside familiar vital sign parameters of ECG and SpO<sub>2</sub>. It is safe, easy to use, economical solution for cardiac output measurement as it has no additional running costs or accessories.



## Visualizing Volumetric Information with Hemodynamics Graph

For quality patient care and comprehensive management of different hemodynamic parameter in septic shock patients.



Nihon Kohden's Hemodynamic Graph provides a more intuitive approach to diagnostic and therapeutic decision making in hemodynamics management. This new tool helps clinicians easily see the direction and trend of hemodynamic change while imaging the Frank Starling curve.



# Improving the Mortality Rate in Sepsis

## Hemodynamics graph for sepsis treatment


Blood Pressure and CVP (Central Venous Pressure) target graph can support Early Goal Directed Therapy (EGDT) according to the guidelines (SSC 2008) for initial resuscitation of severe sepsis and septic shock. The monitor plots the values on the target graph as shown in figure.



CVP/MAP initial resuscitation protocol in severe sepsis and septic shock

## CVP\_ET

Respiration affects intrapleural pressure and CVP value becomes unstable. Nihon Kohden original technology CVP\_ET uses CO<sub>2</sub> signal to calculate the end tidal mean CVP. CVP\_ET provides stable and physiologically - correct value of CVP.




# Effective Alarm System

Life Scope TR easily catches a change in the patient's condition. An effective alarm system and user friendly interface support your daily medical care.

## Interbed monitoring

You can use any bedside monitor to check the vital information and alarm status of another monitor in the network, even if there is no central monitor.

Numeric data for 20 patients or numeric data and 2 waveforms for one patient can be displayed on the Interbed screen.



# High quality monitoring

## Redefining transport

Life Scope TR provides continuous monitoring without reducing the monitoring parameters during transport. When the patient is transported to a new care setting, just detach the input unit and connect it to a different monitor.

The patient information, including trend and waveforms, is transferred to the new bedside monitor to create one seamless patient record.



# External Units



## Ventilators

- Dräger
- Puritan Bennett
- MAQUET
- Hamilton
- GE
- Newport Medical
- Air Liquide
- ResMed
- Metran
- Care Fusion



**Multigas unit**  
GF-210R  
**Multigas/flow unit**  
GF-220R



**Neuro unit**  
AE-918P



**BSM-6000 series**



**BIS® processor**  
QE-910P

## Regional saturation of oxygen (rSO<sub>2</sub>)

- Covidien  
Invos 5100c



## Anesthesia workstation

- Dräger
- GE
- MAQUET
- Heinen + Löwenstein
- Air Liquide

## CCO monitors

- Edwards Lifesciences  
Vigilance, Vigilance II, Vigileo
- Hospira  
Q2™, Q2™ Plus, Q-Vue™  
Pulsion Medical Systems  
PicCO monitor



**NMT Module (TOF)**  
AF-101P

# Major Options



**Input Unit** (3 MULTI connectors)  
AY-663P: for Nihon Kohden (SpO<sub>2</sub>)  
AY-653P: for Nellcor OxiMax® (SpO<sub>2</sub>)  
AY-633P: for Masimo SET® (SpO<sub>2</sub>)

**Input Unit** (1 MULTI connector)  
AY-661P: for Nihon Kohden (SpO<sub>2</sub>)  
AY-651P: for Nellcor OxiMax® (SpO<sub>2</sub>)  
AY-631P: for Masimo SET® (SpO<sub>2</sub>)

**Input Unit** (1 MULTI connector, IBP/CO<sub>2</sub> dedicated type)  
AY-660P: for Nihon Kohden (SpO<sub>2</sub>)

**Memory Unit**  
QM-600P

**Memory Card**  
QM-601P

**Hemodynamics Review Program**  
(including esCCO)  
QP-033P

**esCCO Program**  
QP-034P



**Smart Unit** (4 MULTI connectors)  
AA-674P



**Smart Unit** (2 MULTI connectors)  
AA-672P



**Data Acquisition Unit**  
(no MULTI connectors)  
JA-690PA



**Recorder Module**  
WS-671P



**Remote Controller**  
RY-910PA



**Battery Charger**  
SB-610RK


# Consumables

## CO<sub>2</sub> sensor


Fast response with no lag time  
No warm-up time  
Compact and lightweight

No heater or motor  
No sampling tube  
Durable and easy to use

Mainstream, non-intubated and intubated patients



**cap-ONE**  
ORAL NASAL EXPIRATION



CO<sub>2</sub> sensor kit TG-920P



Airway adapter R804



Oral/nasal adapter V923, with oxygen cannula adapter

## SpO<sub>2</sub> sensor

Water resistant



P225F, finger

Durable




TL-651T3, finger

Durable




P225G, Multisite Type  
Neonate to Adult use

Disposable



P203A, adult finger  
P203B, child finger

# Life Scope Network



This brochure may be revised or replaced by Nihon Kohden at any time without notice.



**NIHON KOHDEN INDIA PVT. LTD.**  
308, Tower-A, Spazedge, Sector-47,  
Sohna Road, Gurgaon 122002 Haryana, India  
Phone +91 124 4931000 (28 Lines) Fax +91 124 4931029  
Email : info\_nki@nkc.co.jp

**NIHON KOHDEN CORPORATION**  
1-31-4 Nishiochiai, Shinjuku-ku, Tokyo 161-8560, Japan  
Phone +81 (3) 5996-8036 Fax +81 (3) 5996-8100  
www.nihonkohden.com

Customer Care Toll Free No.: 1800-103-8182